

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26



Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 2.0 mm; angle-stable; stainless steel</b>		
2.114.06	Cortical screw with conical head thread; Ø 2.0 x 6 mm; fully threaded; self tapping; stainless steel	4 047075 129669
2.114.08	Cortical screw with conical head thread; Ø 2.0 x 8 mm; fully threaded; self tapping; stainless steel	4 047075 129676
2.114.10	Cortical screw with conical head thread; Ø 2.0 x 10 mm; fully threaded; self tapping; stainless steel	4 047075 129683
2.114.12	Cortical screw with conical head thread; Ø 2.0 x 12 mm; fully threaded; self tapping; stainless steel	4 047075 129690
2.114.14	Cortical screw with conical head thread; Ø 2.0 x 14 mm; fully threaded; self tapping; stainless steel	4 047075 129706
2.114.16	Cortical screw with conical head thread; Ø 2.0 x 16 mm; fully threaded; self tapping; stainless steel	4 047075 142170
2.114.18	Cortical screw with conical head thread; Ø 2.0 x 18 mm; fully threaded; self tapping; stainless steel	4 047075 142187
2.114.20	Cortical screw with conical head thread; Ø 2.0 x 20 mm; fully threaded; self tapping; stainless steel	4 047075 129713
2.114.22	Cortical screw with conical head thread; Ø 2.0 x 22 mm; fully threaded; self tapping; stainless steel	4 047075 142194
2.114.24	Cortical screw with conical head thread; Ø 2.0 x 24 mm; fully threaded; self tapping; stainless steel	4 047075 142200
2.114.26	Cortical screw with conical head thread; Ø 2.0 x 26 mm; fully threaded; self tapping; stainless steel	4 047075 142217
<b>Cortical screw Ø 2.7 mm; angle-stable; stainless steel</b>		
2.125.08	Cortical screw with conical head thread; Ø 2.7 x 8 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 129218
2.125.10	Cortical screw with conical head thread; Ø 2.7 x 10 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124244
2.125.12	Cortical screw with conical head thread; Ø 2.7 x 12 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124251
2.125.14	Cortical screw with conical head thread; Ø 2.7 x 14 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124268
2.125.16	Cortical screw with conical head thread; Ø 2.7 x 16 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124275

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26



Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 2.7 mm; angle-stable; stainless steel</b>		
2.125.18	Cortical screw with conical head thread; Ø 2.7 x 18 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124282
2.125.20	Cortical screw with conical head thread; Ø 2.7 x 20 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124299
2.125.22	Cortical screw with conical head thread; Ø 2.7 x 22 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124305
2.125.24	Cortical screw with conical head thread; Ø 2.7 x 24 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124312
2.125.26	Cortical screw with conical head thread; Ø 2.7 x 26 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124329
2.125.28	Cortical screw with conical head thread; Ø 2.7 x 28 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124336
2.125.30	Cortical screw with conical head thread; Ø 2.7 x 30 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 124343
2.125.32	Cortical screw with conical head thread; Ø 2.7 x 32 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161683
2.125.34	Cortical screw with conical head thread; Ø 2.7 x 34 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161690
2.125.36	Cortical screw with conical head thread; Ø 2.7 x 36 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161706
2.125.38	Cortical screw with conical head thread; Ø 2.7 x 38 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161713
2.125.40	Cortical screw with conical head thread; Ø 2.7 x 40 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161720
2.125.42	Cortical screw with conical head thread; Ø 2.7 x 42 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161737
2.125.44	Cortical screw with conical head thread; Ø 2.7 x 44 mm; fully threaded; self tapping; head 5.0 mm; stainless steel	4 047075 161744

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26




















Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 3.5 mm; angle-stable; stainless steel</b>		
2.133.10	Cortical screw with conical head thread; Ø 3.5 x 10 mm; fully threaded; self tapping; stainless steel	
2.133.12	Cortical screw with conical head thread; Ø 3.5 x 12 mm; fully threaded; self tapping; stainless steel	
2.133.14	Cortical screw with conical head thread; Ø 3.5 x 14 mm; fully threaded; self tapping; stainless steel	
2.133.16	Cortical screw with conical head thread; Ø 3.5 x 16 mm; fully threaded; self tapping; stainless steel	
2.133.18	Cortical screw with conical head thread; Ø 3.5 x 18 mm; fully threaded; self tapping; stainless steel	
2.133.20	Cortical screw with conical head thread; Ø 3.5 x 20 mm; fully threaded; self tapping; stainless steel	
2.133.22	Cortical screw with conical head thread; Ø 3.5 x 22 mm; fully threaded; self tapping; stainless steel	
2.133.24	Cortical screw with conical head thread; Ø 3.5 x 24 mm; fully threaded; self tapping; stainless steel	
2.133.26	Cortical screw with conical head thread; Ø 3.5 x 26 mm; fully threaded; self tapping; stainless steel	
2.133.28	Cortical screw with conical head thread; Ø 3.5 x 28 mm; fully threaded; self tapping; stainless steel	
2.133.30	Cortical screw with conical head thread; Ø 3.5 x 30 mm; fully threaded; self tapping; stainless steel	
2.133.32	Cortical screw with conical head thread; Ø 3.5 x 32 mm; fully threaded; self tapping; stainless steel	
2.133.34	Cortical screw with conical head thread; Ø 3.5 x 34 mm; fully threaded; self tapping; stainless steel	
2.133.35	Cortical screw with conical head thread; Ø 3.5 x 35 mm; fully threaded; self tapping; stainless steel	
2.133.36	Cortical screw with conical head thread; Ø 3.5 x 36 mm; fully threaded; self tapping; stainless steel	
2.133.38	Cortical screw with conical head thread; Ø 3.5 x 38 mm; fully threaded; self tapping; stainless steel	
2.133.40	Cortical screw with conical head thread; Ø 3.5 x 40 mm; fully threaded; self tapping; stainless steel	

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26



Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 3.5 mm; angle-stable; stainless steel</b>		
2.133.42	Cortical screw with conical head thread; Ø 3.5 x 42 mm; fully threaded; self tapping; stainless steel	4 047075 069446
2.133.44	Cortical screw with conical head thread; Ø 3.5 x 44 mm; fully threaded; self tapping; stainless steel	4 047075 069347
2.133.45	Cortical screw with conical head thread; Ø 3.5 x 45 mm; fully threaded; self tapping; stainless steel	4 047075 113996
2.133.46	Cortical screw with conical head thread; Ø 3.5 x 46 mm; fully threaded; self tapping; stainless steel	4 047075 069354
2.133.48	Cortical screw with conical head thread; Ø 3.5 x 48 mm; fully threaded; self tapping; stainless steel	4 047075 069255
2.133.50	Cortical screw with conical head thread; Ø 3.5 x 50 mm; fully threaded; self tapping; stainless steel	4 047075 069262
2.133.55	Cortical screw with conical head thread; Ø 3.5 x 55 mm; fully threaded; self tapping; stainless steel	4 047075 113972
2.133.60	Cortical screw with conical head thread; Ø 3.5 x 60 mm; fully threaded; self tapping; stainless steel	4 047075 069378
2.133.65	Cortical screw with conical head thread; Ø 3.5 x 65 mm; fully threaded; self tapping; stainless steel	4 047075 069293
<b>Cortical screw Ø 4.5 mm; angle-stable; stainless steel</b>		
2.545.10	Cortical screw with conical head thread; Ø 4.5 x 10 mm; fully threaded; self tapping; stainless steel	4 047075 080748
2.545.12	Cortical screw with conical head thread; Ø 4.5 x 12 mm; fully threaded; self tapping; stainless steel	4 047075 080632
2.545.14	Cortical screw with conical head thread; Ø 4.5 x 14 mm; fully threaded; self tapping; stainless steel	4 047075 080649
2.545.16	Cortical screw with conical head thread; Ø 4.5 x 16 mm; fully threaded; self tapping; stainless steel	4 047075 080779
2.545.18	Cortical screw with conical head thread; Ø 4.5 x 18 mm; fully threaded; self tapping; stainless steel	4 047075 080878
2.545.20	Cortical screw with conical head thread; Ø 4.5 x 20 mm; fully threaded; self tapping; stainless steel	4 047075 080885
2.545.22	Cortical screw with conical head thread; Ø 4.5 x 22 mm; fully threaded; self tapping; stainless steel	4 047075 080755

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26



Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 4.5 mm; angle-stable; stainless steel</b>		
2.545.24	Cortical screw with conical head thread; Ø 4.5 x 24 mm; fully threaded; self tapping; stainless steel	4 047075 080762
2.545.26	Cortical screw with conical head thread; Ø 4.5 x 26 mm; fully threaded; self tapping; stainless steel	4 047075 080915
2.545.28	Cortical screw with conical head thread; Ø 4.5 x 28 mm; fully threaded; self tapping; stainless steel	4 047075 080922
2.545.30	Cortical screw with conical head thread; Ø 4.5 x 30 mm; fully threaded; self tapping; stainless steel	4 047075 080793
2.545.32	Cortical screw with conical head thread; Ø 4.5 x 32 mm; fully threaded; self tapping; stainless steel	4 047075 080809
2.545.34	Cortical screw with conical head thread; Ø 4.5 x 34 mm; fully threaded; self tapping; stainless steel	4 047075 080892
2.545.36	Cortical screw with conical head thread; Ø 4.5 x 36 mm; fully threaded; self tapping; stainless steel	4 047075 080908
2.545.38	Cortical screw with conical head thread; Ø 4.5 x 38 mm; fully threaded; self tapping; stainless steel	4 047075 081028
2.545.40	Cortical screw with conical head thread; Ø 4.5 x 40 mm; fully threaded; self tapping; stainless steel	4 047075 080939
2.545.42	Cortical screw with conical head thread; Ø 4.5 x 42 mm; fully threaded; self tapping; stainless steel	4 047075 080946
2.545.44	Cortical screw with conical head thread; Ø 4.5 x 44 mm; fully threaded; self tapping; stainless steel	4 047075 080816
2.545.45	Cortical screw with conical head thread; Ø 4.5 x 45 mm; fully threaded; self tapping; stainless steel	4 047075 080823
2.545.46	Cortical screw with conical head thread; Ø 4.5 x 46 mm; fully threaded; self tapping; stainless steel	4 047075 080977
2.545.48	Cortical screw with conical head thread; Ø 4.5 x 48 mm; fully threaded; self tapping; stainless steel	4 047075 081035
2.545.50	Cortical screw with conical head thread; Ø 4.5 x 50 mm; fully threaded; self tapping; stainless steel	4 047075 081042
2.545.52	Cortical screw with conical head thread; Ø 4.5 x 52 mm; fully threaded; self tapping; stainless steel	4 047075 080953
2.545.54	Cortical screw with conical head thread; Ø 4.5 x 54 mm; fully threaded; self tapping; stainless steel	4 047075 080960

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26



Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 4.5 mm; angle-stable; stainless steel</b>		
2.545.55	Cortical screw with conical head thread; Ø 4.5 x 55 mm; fully threaded; self tapping; stainless steel	4 047075 081073
2.545.56	Cortical screw with conical head thread; Ø 4.5 x 56 mm; fully threaded; self tapping; stainless steel	4 047075 081080
2.545.58	Cortical screw with conical head thread; Ø 4.5 x 58 mm; fully threaded; self tapping; stainless steel	4 047075 080984
2.545.60	Cortical screw with conical head thread; Ø 4.5 x 60 mm; fully threaded; self tapping; stainless steel	4 047075 080991
2.545.62	Cortical screw with conical head thread; Ø 4.5 x 62 mm; fully threaded; self tapping; stainless steel	4 047075 081059
2.545.64	Cortical screw with conical head thread; Ø 4.5 x 64 mm; fully threaded; self tapping; stainless steel	4 047075 081066
2.545.65	Cortical screw with conical head thread; Ø 4.5 x 65 mm; fully threaded; self tapping; stainless steel	4 047075 033065
2.545.66	Cortical screw with conical head thread; Ø 4.5 x 66 mm; fully threaded; self tapping; stainless steel	4 047075 081189
2.545.68	Cortical screw with conical head thread; Ø 4.5 x 68 mm; fully threaded; self tapping; stainless steel	4 047075 081097
2.545.70	Cortical screw with conical head thread; Ø 4.5 x 70 mm; fully threaded; self tapping; stainless steel	4 047075 081103
2.545.75	Cortical screw with conical head thread; Ø 4.5 x 75 mm; fully threaded; self tapping; stainless steel	4 047075 069972
2.545.80	Cortical screw with conical head thread; Ø 4.5 x 80 mm; fully threaded; self tapping; stainless steel	4 047075 094424
2.545.85	Cortical screw with conical head thread; Ø 4.5 x 85 mm; fully threaded; self tapping; stainless steel	4 047075 094431

# Cortical screw

## angle-stable; stainless steel



Item no. 2.114.14



Item no. 2.125.26



Item no. 2.133.20



Item no. 2.545.24



Item no. 2.549.40

Item no.	Item description	EAN
<b>Cortical screw Ø 5.5 mm; angle-stable; stainless steel</b>		
2.549.20	Cortical screw with conical head thread; Ø 5.5 x 20 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 146970
2.549.35	Cortical screw with conical head thread; Ø 5.5 x 35 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 146987
2.549.40	Cortical screw with conical head thread; Ø 5.5 x 40 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139569
2.549.45	Cortical screw with conical head thread; Ø 5.5 x 45 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139576
2.549.50	Cortical screw with conical head thread; Ø 5.5 x 50 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139583
2.549.55	Cortical screw with conical head thread; Ø 5.5 x 55 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139590
2.549.60	Cortical screw with conical head thread; Ø 5.5 x 60 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139606
2.549.65	Cortical screw with conical head thread; Ø 5.5 x 65 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139613
2.549.70	Cortical screw with conical head thread; Ø 5.5 x 70 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139620
2.549.75	Cortical screw with conical head thread; Ø 5.5 x 75 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139637
2.549.80	Cortical screw with conical head thread; Ø 5.5 x 80 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139644
2.549.85	Cortical screw with conical head thread; Ø 5.5 x 85 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139651
2.549.90	Cortical screw with conical head thread; Ø 5.5 x 90 mm; fully threaded; self tapping; stainless steel	4 <sup>0</sup> 047075 <sup>1</sup> 139668