


## All advantages at a glance

- Short and long nails; various diameters (10 mm - 13 mm)
- Various CCD angles (125° and 130°)
- Colour-coded instruments for easier handling
- Special supporting screw design with controlled compression path of 15 mm, providing in a high degree of primary stability
- Radiolucent, lightweight, carbon fibre-reinforced plastic guide frame with modular technology



## Contact us

If we have piqued your interest in the KFN-Königsee Femoral Nail, please contact us or your sales representative.



### Königsee Implantate GmbH

OT Aschau | Am Sand 4  
07426 Allendorf  
Germany

export@koenigsee-implantate.de  
+ 49 36738 498-560

ALL PRODUCTS ALSO AVAILABLE  
UNDER THE MDR IN NON-STERILE  
AND STERILE.



## KFN-KÖNIGSEE FEMORAL NAIL

### Proximal femoral nail



[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)


Copyright © Königsee Implantate GmbH  
SiZi | Edition 03 | 2024-06



 **Königsee**  
Implantate

## Target group

The KFN is intended for adult patients. It may not be used for children or young people who have not finished growing.



## Properties

### Product properties - short nails:

- Length: 210 mm
- Cannulated
- Titanium alloy
- anodised surface

### Colour coding of nail diameters:

- Ø 10 mm: anodised purple
- Ø 11 mm: anodised green-blue
- Ø 12 mm: anodised blue
- Ø 13 mm: anodised brown

### Product variants:

- Item no. 7.1KF.XX/210
- CCD angle: 130°
- Item no. 7.3KF.XX/210
- CCD angle: 125°


### Product properties - long nails:


- Ø 10.5 mm
- Cannulated
- Titanium alloy
- anodised golden yellow

### Product variants:

- Item no. 7.1KF10,5L/XXX (left femur)
- Item no. 7.1KF10,5R/XXX (right femur)
- Lengths: 340 mm, 420 mm
- CCD angle: 130°
- Item no. 7.3KF10,5L/XXX (left femur)
- Item no. 7.3KF10,5R/XXX (right femur)
- Lengths: 340 mm, 380 mm
- CCD angle: 125°

## Set Königsee Femoral Nail



Set implants; Set no. 19.344.00  Internal hexagon  
Container implants; Item no. 19.345.00




Item no.  
7.1KF.10/210


Item no.  
7.1KF.11/210

Item no.  
7.1KF.12/210

Container 1 instruments; Item no. 19.347.00  
Set 1 instruments; Set no. 19.344.00



Container 2 instruments; Item no. 19.348.00  
Set 2 instruments; Set no. 19.344.00



Container 3 instruments; Item no. 19.349.00  
Set 3 instruments; Set no. 19.344.00

