

# KFN-Königsee Femoral Nail

Titanium

not sterile packed / sterile packed



Item no. 7.1KF.10/210  
Item no. sterile S7.1KF.10/210



Item no. 7.1KF.11/210  
Item no. sterile S7.1KF.11/210



Item no. 7.1KF.12/210  
Item no. sterile S7.1KF.12/210

Item no.	EAN	Item description	Item no. sterile	EAN sterile
----------	-----	------------------	------------------	-------------

## Implants - Nails

### Königsee Femoral Nail 125°; titanium

7.3KF.10/210		Königsee Femoral Nail 125°; Ø 10 x 210 mm; cannulated; inside Ø 4.7 mm	S7.3KF.10/210	
7.3KF.11/210		Königsee Femoral Nail 125°; Ø 11 x 210 mm; cannulated; inside Ø 4.7 mm	S7.3KF.11/210	
7.3KF.12/210		Königsee Femoral Nail 125°; Ø 12 x 210 mm; cannulated; inside Ø 4.7 mm	S7.3KF.12/210	
7.3KF.13/210		Königsee Femoral Nail 125°; Ø 13 x 210 mm; cannulated; inside Ø 4.7 mm	S7.3KF.13/210	
7.3KF10,5L/340		Königsee Femoral Nail 125°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; left	S7.3KF.340L	
7.3KF10,5R/340		Königsee Femoral Nail 125°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; right	S7.3KF.340R	
7.3KF10,5L/380		Königsee Femoral Nail 125°; Ø 10.5 x 380 mm; cannulated; inside Ø 4.7 mm; left	S7.3KF.380L	
7.3KF10,5R/380		Königsee Femoral Nail 125°; Ø 10.5 x 380 mm; cannulated; inside Ø 4.7 mm; right	S7.3KF.380R	

### Königsee Femoral Nail 130°; titanium

7.1KF.10/210		Königsee Femoral Nail 130°; Ø 10 x 210 mm; cannulated; inside Ø 4.7 mm	S7.1KF.10/210	
7.1KF.11/210		Königsee Femoral Nail 130°; Ø 11 x 210 mm; cannulated; inside Ø 4.7 mm	S7.1KF.11/210	
7.1KF.12/210		Königsee Femoral Nail 130°; Ø 12 x 210 mm; cannulated; inside Ø 4.7 mm	S7.1KF.12/210	
7.1KF.13/210		Königsee Femoral Nail 130°; Ø 13 x 210 mm; cannulated; inside Ø 4.7 mm	S7.1KF.13/210	
7.1KF10,5L/340		Königsee Femoral Nail 130°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; left	S7.1KF.340L	
7.1KF10,5R/340		Königsee Femoral Nail 130°; Ø 10.5 x 340 mm; cannulated; inside Ø 4.7 mm; right	S7.1KF.340R	
7.1KF10,5L/420		Königsee Femoral Nail 130°; Ø 10.5 x 420 mm; cannulated; inside Ø 4.7 mm; left	S7.1KF.420L	
7.1KF10,5R/420		Königsee Femoral Nail 130°; Ø 10.5 x 420 mm; cannulated; inside Ø 4.7 mm; right	S7.1KF.420R	



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany

Fon +49 36738 498-560  
Fax +49 36738 498-569

Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

For further information please visit:  
[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)

StZi | Edition 03 | 2024-06  
All pictured products are CE compliant.

# KFN-Königsee Femoral Nail

Titanium

not sterile packed / sterile packed



Item no. 3.157.36  
Item no. sterile S3.157.36

Item no.	EAN	Item description	Item no. sterile	EAN sterile
<b>Implants - Screws; titanium alloy; TiOB<sup>®</sup>-coated</b>				
3.157.25		interlocking screw for KFN; Ø 5.5 x 25 mm; self tapping	S3.157.25	
3.157.26		interlocking screw for KFN; Ø 5.5 x 26 mm; self tapping	S3.157.26	
3.157.28		interlocking screw for KFN; Ø 5.5 x 28 mm; self tapping	S3.157.28	
3.157.30		interlocking screw for KFN; Ø 5.5 x 30 mm; self tapping	S3.157.30	
3.157.32		interlocking screw for KFN; Ø 5.5 x 32 mm; self tapping	S3.157.32	
3.157.34		interlocking screw for KFN; Ø 5.5 x 34 mm; self tapping	S3.157.34	
3.157.35		interlocking screw for KFN; Ø 5.5 x 35 mm; self tapping	S3.157.35	
3.157.36		interlocking screw for KFN; Ø 5.5 x 36 mm; self tapping	S3.157.36	
3.157.38		interlocking screw for KFN; Ø 5.5 x 38 mm; self tapping	S3.157.38	
3.157.40		interlocking screw for KFN; Ø 5.5 x 40 mm; self tapping	S3.157.40	
3.157.42		interlocking screw for KFN; Ø 5.5 x 42 mm; self tapping	S3.157.42	
3.157.44		interlocking screw for KFN; Ø 5.5 x 44 mm; self tapping	S3.157.44	
3.157.45		interlocking screw for KFN; Ø 5.5 x 45 mm; self tapping	S3.157.45	
3.157.50		interlocking screw for KFN; Ø 5.5 x 50 mm; self tapping	S3.157.50	
3.157.55		interlocking screw for KFN; Ø 5.5 x 55 mm; self tapping	S3.157.55	
3.157.60		interlocking screw for KFN; Ø 5.5 x 60 mm; self tapping	S3.157.60	
3.157.65		interlocking screw for KFN; Ø 5.5 x 65 mm; self tapping	S3.157.65	

# KFN-Königsee Femoral Nail

## Titanium

not sterile packed / sterile packed



Item no. 3.428.100  
Item no. sterile S3.428.100



Item no. 7.311.12  
Item no. sterile S7.311.12



Item no. 7.311.30  
Item no. sterile S7.311.30



Item no. 3.310.52  
Item no. sterile S3.310.52



Item no. 7.311.01  
Item no. sterile S7.311.01

Item no.	EAN	Item description	Item no. sterile	EAN sterile
----------	-----	------------------	------------------	-------------

### Implants - Screws; titanium alloy; TiOB®-coated

3.428.65		anti-rotation screw for KFN; Ø 6.0 x 65 mm	S3.428.65	
3.428.70		anti-rotation screw for KFN; Ø 6.0 x 70 mm	S3.428.70	
3.428.75		anti-rotation screw for KFN; Ø 6.0 x 75 mm	S3.428.75	
3.428.80		anti-rotation screw for KFN; Ø 6.0 x 80 mm	S3.428.80	
3.428.85		anti-rotation screw for KFN; Ø 6.0 x 85 mm	S3.428.85	
3.428.90		anti-rotation screw for KFN; Ø 6.0 x 90 mm	S3.428.90	
3.428.95		anti-rotation screw for KFN; Ø 6.0 x 95 mm	S3.428.95	
3.428.100		anti-rotation screw for KFN; Ø 6.0 x 100 mm	S3.428.100	
3.428.105		anti-rotation screw for KFN; Ø 6.0 x 105 mm	S3.428.105	
3.428.110		anti-rotation screw for KFN; Ø 6.0 x 110 mm	S3.428.110	
3.428.115		anti-rotation screw for KFN; Ø 6.0 x 115 mm	S3.428.115	
3.428.120		anti-rotation screw for KFN; Ø 6.0 x 120 mm	S3.428.120	

### Implants - Screws

7.311.12		locking screw without arresting bolt for KFN; length 28 mm; M 12; titanium alloy; TiOB®-coated	S7.311.12	
7.311.30		locking screw with arresting bolt for KFN; length 30 mm; M 12; titanium alloy; TiOB®-coated	S7.311.30	
7.311.01		compression screw for DHS; M 4 x 27 mm; length 35 mm; SW 3.5; titanium alloy	S7.311.01	
3.310.52		supporting screw sleeve for KFN; Ø 12 x 52 mm; titanium alloy; TiOB®-coated	S3.310.52	

### Implants - Screws; titanium alloy; TiOB®-coated

7.313.70		supporting screw for KFN; Ø 10 x 70 mm	S7.313.70	
7.313.75		supporting screw for KFN; Ø 10 x 75 mm	S7.313.75	























# KFN-Königsee Femoral Nail

Titanium

not sterile packed / sterile packed



Item no. 7.313.120  
Item no. sterile S7.313.120

Item no.	EAN	Item description	Item no. sterile	EAN sterile
<b>Implants - Screws; titanium alloy; TiOB®-coated</b>				
7.313.80		supporting screw for KFN; Ø 10 x 80 mm	S7.313.80	
7.313.85		supporting screw for KFN; Ø 10 x 85 mm	S7.313.85	
7.313.90		supporting screw for KFN; Ø 10 x 90 mm	S7.313.90	
7.313.95		supporting screw for KFN; Ø 10 x 95 mm	S7.313.95	
7.313.100		supporting screw for KFN; Ø 10 x 100 mm	S7.313.100	
7.313.105		supporting screw for KFN; Ø 10 x 105 mm	S7.313.105	
7.313.110		supporting screw for KFN; Ø 10 x 110 mm	S7.313.110	
7.313.115		supporting screw for KFN; Ø 10 x 115 mm	S7.313.115	
7.313.120		supporting screw for KFN; Ø 10 x 120 mm	S7.313.120	
7.313.125		supporting screw for KFN; Ø 10 x 125 mm	S7.313.125	
7.313.130		supporting screw for KFN; Ø 10 x 130 mm	S7.313.130	

## Information

TiOB® is a surface treatment that brings optimized properties.



Königsee Implantate GmbH  
OT Aschau | Am Sand 4  
07426 Allendorf | Germany

Fon +49 36738 498-560  
Fax +49 36738 498-569

Mail [export@koenigsee-implantate.de](mailto:export@koenigsee-implantate.de)

For further information please visit:  
[www.koenigsee-implantate.de](http://www.koenigsee-implantate.de)

StZi | Edition 03 | 2024-06  
All pictured products are CE compliant.