

All advantages at a glance

- Percutaneous implantation for **optimal protection of soft tissue**
- **Stable cortical implant anchoring** regardless of drilling duration
- Fracture compression
- Early functional follow-up treatment is possible after implanting **just one compression wire**, provided fracture healing and bone quality are sufficient
- No immobilisation required after surgery
- Material is easy to remove

Contact details

If you are interested in **PERCUTANEOUS COMPRESSION WIRE** please feel free to contact us or your local sales representatives.



Königsee Implantate GmbH


OT Aschau | Am Sand 4
07426 Allendorf
Germany

export@koenigsee-implantate.de
+ 49 36738 498-560

All implants are available on request as non-sterile packed or sterile packed.

PERCUTANEOUS COMPRESSION WIRE

for the treatment of
finger shaft fractures



www.koenigsee-implantate.de

Copyright © Königsee Implantate GmbH/HeBo | Edition 01 | 2019-06

 **Königsee**
Implantate

Application possibilities

- Simple transverse and short oblique fractures of the proximal and intermediate phalanges

Target group

- The percutaneous compression wire is intended for adolescent and adult patients.
- There are no restrictions with regard to selected ethnic groups. The patient clientèle is only limited by anatomical parameters.

Features

- Implant steel with electropolished surface
- \varnothing 1.2 mm / 1.8 mm
- Length 150 mm
- Distance between threads 5 mm, 10 mm or 15 mm

