

All advantages at a glance

- Anatomically formed, variable angle stable implant
- Uniform plate versions for left and right
- Good mouldability of the implant
- Bilateral fit minimises soft tissue irritation
- Two small arms to facilitate support of the distal pole of the patella
- The angle stable bilateral plate osteosynthesis offers a higher level of stability and rigidity compared to conventional treatment with tension band wiring or cannulated screws.
- Reduction of the risk of secondary dislocation
- Primary reduction result is maintained through the fixed angle principle of fixateur interne

Contact details

If you are interested in the **patella plate 3.5 vas** please feel free to contact us or your local sales representatives.



Königsee Implantate GmbH

OT Aschau | Am Sand 4
07426 Allendorf
Germany

export@koenigsee-implantate.de
+ 49 36738 498 - 560

All implants are available on request as non-sterile packed or sterile packed.

PATELLA PLATE

Variable angle stable patella plate 3.5



www.koenigsee-implantate.de

Target group

Treatment with a variable angle stable patella plate 3.5 is carried out on adolescent and adult patients.

Features

Variable angle stable patella plate 3.5

Pure titanium

Blue anodised

11 plate holes


Anatomically adapted plate shape

Uniform plate versions for left and right

Variable angle stable holes for angle stable cortical screws \varnothing 3.5 mm; head \varnothing 5.0 mm

Eyelets for temporary fixation

Set Patella plate 3.5 vas



OP set patella plate 3.5 vas with instruments and screws

 Internal hexagon

